

SAFE RMR TESTING

DURING THE COVID-19 PANDEMIC

1 PROTECT PATIENTS USE A KORR METABREATHER

KORR MetaBreathers are **DISPOSABLE** and **STERILE**, protecting patients with a **ONE-WAY VALVE** that draws fresh air from the room with each breath.

2 PROTECT AGAINST CROSS-CONTAMINATION USE AN ANTI-VIRAL FILTER

PROTECT YOUR DEVICE from infected patients with an anti-viral filter. They're FDA cleared to be **99.9% EFFECTIVE** and usually **LAST FOR 100 TESTS**.

3 PROTECT HEALTHCARE WORKERS PROPER PROTOCOL & PPE

WEAR MASKS AND GLOVES when disposing of used MetaBreather tubing. **DON'T TEST PATIENTS WHO ARE SYMPTOMATIC** or who test positive for COVID-19.

WHY IS TESTING SAFE?

THE DISPOSABLE MOUTHPIECE TOUCHING YOUR FACE IS STERILE

ONLY CLEAN AIR ENTERS THE ONE-WAY VALVE

NO CROSS-CONTAMINATION BECAUSE THE AIR MOVES ONE WAY

